




Water Still

Full Safe
High performance



Technical characteristic

- Durable, reliable and Low cost.
- High purity distilled water - conductivity approx 2,5 $\mu\text{S}/\text{cm}$.
- External tank level sensor.
- Automatic water switching and power cut-off.
- Automatic thermostatic cut-off safety system.
- Heating element made of high quality stainless steel.
- All material contact water made of stainless steel 304L.
- Easily accessed evaporator & storage tank for effortless cleaning and maintenance.
- The unit is suitable for both bench and wall mounting.
- Available Pro version with advance diagnostic system detect heating element failed

	Capacity	Consumption	Power	Power Requirement	Weight	Dimension(mm)
WS-02L	2L/H	25L/H	1700w	AC 220V	15.8kg	210x250x410
WS-04L	4L/H	40L/H	3000w	AC 220V	20.2kg	280x250x520
WS-08L	8L/H	70L/H	6000w	3Ph AC 380V	32.5kg	320*450*620
WS-12L	12L/H	193L/H	9000w	3Ph AC 380V	44.2kg	450*650*700

2 year warranty (heating element not include)



Contact

Facebook page : facebook.com/rlabtn

Website : www.rlab.tn

Tunisia office: 

Phone: +216 73 534 458

Mobile: +216 93 413 441

Technical support : tech@rlab.tn

Sales department: sales@rlab.tn

Information service: info@rlab.tn