


Water Bath

High performance and reliability.
Simplicity. Efficacy and saving cost.



Technical characteristic

- PID controller with Auto-tuning provides accurate temperature regulation.
- Case mild steel with anti-corrosion treatment and epoxy paint.
- Temperature range (Ambient T + 5°C - 70°C) without lid.
- Temperature range (Ambient T + 5°C - 100°C) with lid.
- Thermal isolation Rock-wool 30mm high density.
- Temperature Sensor PT100 class A.
- Temperature uniformity: $\pm 0.2^{\circ}\text{C}$ @ 37.0°C.
- Temperature Stability: $\pm 1.0^{\circ}\text{C}$.
- Bath made from Stainless Steel.
- Over temperature protection.

	WB-11	WB-17	WB-34
Reservoir Capacity (L)	11	17	34
Bath Dimensions (L x W x D mm)	170 x 320 x 200	260 x 320 x 200	530 x 320 x 200
Overall Dimensions (L x W x D mm)	270 x 400 x 230	320 x 400 x 230	600 x 420 x 230
Power (W)	500	700	1200
Power Requirement	AC 220V	AC 220V	AC 220V
Weight (Kg)	10.8	18.6	27.3

2 year warranty



Contact

Facebook page : facebook.com/rlabtn

Website : www.rlab.tn

Tunisia office: 

Phone: +216 73 534 458

Mobile: +216 93 413 441

Technical support : tech@rlab.tn

Sales department: sales@rlab.tn

Information service: info@rlab.tn